



## Prepurification of adipokinetic hormones

Adipokinetic hormones (AKHs) are insect anti-stress hormones that maintain the biochemical and physiological homeostasis of the insect body (Kodrík, 2008). AKHs are octa, nona- or decapeptides with both termini blocked: the N-terminus by a pyroglutamate residue and the C-terminus by an amide. Typically, a specific antibody and ELISA method are used for their quantification in the insect central nervous system and haemolymph. For the latter, prior to executing the ELISA test, several pre-purification steps are required; they also involve a solid phase extraction cartridge.

In the test, the AKH from the firebug *Pyrrhocoris apterus* known as Pyrap-AKH was employed. Its structure is: pGlu-Leu-Asn-Phe-Thr-Pro-Asn-Trp-NH<sub>2</sub> (Kodrík et al., 2000).



*Sample Separation – 160 pmol-Pyrap-AKH*



*Control Sample – 160 pmol-Pyrap-AKH*

### References

- Kodrík D. (2008) Adipokinetic hormone functions that are not associated with insect flight. *Physiol. Entomol.* 33: 171-180.
- Kodrík D., Socha R., Šimek P., Zemek R. and Goldsworthy G.J. (2000) A new member of the AKH/RPCH family that stimulates locomotory activity in the firebug, *Pyrrhocoris apterus* (Heteroptera), *Insect Biochem. Mol. Biol.* 30: 489-498.



## Prepurification of adipokinetic hormones

### MSPE Method

<b>MSP column</b>	Micro Spin SpeExtra™ column MSPE C18P 0.22 µm Nylon membrane, 15 mg, 0.7 mL
<b>Solution A</b>	0.11% trifluoroacetic acid (TFA) in water
<b>Solution B</b>	0.1% TFA in 60% acetonitrile
<b>Centrifugation</b>	2000 rev/min
<b>MSPE steps</b>	<ol style="list-style-type: none"><li><b>1. Solution B, 0.6 mL</b></li><li><b>2. Solution A, 0.6 mL</b></li><li><b>3. Pyrap-AKH 160 pmol in Solution A, 0.6 mL</b></li><li><b>4. The eluate applied again on the cartridge, 0.6 mL</b></li><li><b>5. Solution A, 0.6 mL</b></li><li><b>6. Elution with 0.3 mL Solution B and 0.3 mL 100% acetonitrile</b></li></ol>

### HPLC method

<b>Column</b>	Chromolith Performance RP-18e column 150 × 4.6 mm (Merck)
<b>Mobile phase</b>	A = 0.11% TFA in water, B = 0.1% TFA in 60% acetonitrile
<b>Gradient</b>	0–2 min 30 % B, 2–12 min 30–90 % B, 12–15 min 90 % B
<b>Flow rate</b>	1.5 mL/min
<b>Detection</b>	Fluorescence, wavelength Ex 280 and Em 348 nm (Waters model 2475)

This application was developed  
by the Biology centre CAS, České  
Budějovice, The Czech Republic

