



## Astra® C18-Hybrid ChromLine

# ESTROGENS ANALYSIS WITH ASTRA® C18-HYBRID COLUMN

### ABSTRACT


This application describes a quick and robust LC-MS/MS method for analysis of estrone (E1 – shown in the upper chromatograms) and estradiol (E2 – shown in the lower chromatograms), two naturally occurring estrogens in the human body. The ASTRA® C18-Hybrid column provides perfect selectivity and peak symmetry for target compounds throughout a long injection sequence.

### INTRODUCTION


Analysis of endogenous estrogens is challenging due to the low physiological concentrations and weak UV absorption. Multiple metabolites with similar chemical structures and close molecular weights as well as co-eluting substances from biological fluids may cause significant ion suppression in mass spectrometry. Ionization of estrogens is inefficient, introduction of derivatization agents may enhance ionization, but non-specific product ions can be generated. Therefore, the method generating derivatives with product ions specific to individual estrogens is highly demanded.

Astra® C18-Hybrid  
ChromLine


## Estrogens Analysis with ASTRA® C18-Hybrid column




Analysis of derivatized hormones E1 and E2 in standard mixture and plasma sample on ASTRA® C18-Hybrid column



Analysis of plasma sample on different lots of ASTRA® C18-Hybrid column

**Astra® C18-Hybrid**  
ChromLine**Estrogens Analysis with ASTRA® C18-Hybrid column**



*Analysis of plasma sample on new column ASTRA® C18-Hybrid and after 1000 injections*



*Analysis of plasma sample on column ASTRA® C18-Hybrid and competitive column*

**Astra® C18-Hybrid**  
ChromLine**Estrogens Analysis with ASTRA® C18-Hybrid column**

<b>Column</b>	ASTRA® C18-Hybrid, 5 µm
<b>Dimensions</b>	50 mm × 4.6 mm
<b>Part number</b>	AST-6016-LG46
<b>Mobile phase</b>	70 % MeOH + 10 % 50mM (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub> + 20 % UPW Detailed information about gradient table is proprietary
<b>Flow rate</b>	1 mL/min
<b>Temperature</b>	50 °C
<b>Injection volume</b>	60 µL
<b>LC-MS conditions</b>	Proprietary
<b>Analytes</b>	<b>1. Estrone (E1), CAS: 53-16-7</b> <b>2. Estradiol (E2), CAS: 50-28-2</b>

*Estrone (E1)**Estradiol (E2)*



The premium product brand of Chromservis s.r.o.



Jakobiho 327/3 | 109 00 Praha 10 – Petrovice | Czech Republic | [www.chromline.eu](http://www.chromline.eu)

---

**ASTRA®** is a trademark by Chromservis s.r.o. Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners. We disclaim proprietary interest in the marks and names of others.