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USER MANUAL

SENSORS CHS SERIES FOR CONDUCTIVITY MEASUREMENT

Sensor overview

This instruction is intended for use with following sensors:
CHS CondiGo, CHS CondiLab

First use

Check that the probe measurement range and cell constant are appropriate for your application.

Depending on the instrument type, it may be necessary to enter the cell constant of the probe. See the table below.

Rinse the sensor with distilled water before first use. There is usually no need for more preparations. However, if measurement shows a larger deviation, immerse the probe in distilled water for a few minutes.

If the sensor body is made from glass, be careful when handling it.

Connect the sensor to the device. When connecting / disconnecting the cable, always handle only the connector, not the cable! It may tear.

Measurement

Perform calibration (see below) before measuring samples. To avoid sample contamination, the sensor should be rinsed in distilled or deionized water when switching between samples

Bubbles can significantly affect measurement. If a bubble appears in the cell, remove the probe from the sample and re-immersing it.

Always immerse the entire measurement cell in the sample! The minimum probe immersion depth is in the table below.

Sensor calibration

The instrument-probe assembly needs to be periodically calibrated to reach required measurement accuracy. The calibration interval depends on the quality standards of the workplace. For calibration, follow the instructions of the conductometer. The general procedure is: rinse the probe with distilled water, then rinse it with the calibration solution, and after that, immerse the probe in the calibration solution. Use standards from lower to higher values to reduce the possibility of contamination of low-conductivity standards.

For calibration solutions below 84 µS/cm, follow the instructions for the solution, usually calibrate directly in the storage bottle.

Sensor storage

Rinse the sensor with clean water after use. The probe can be stored both in distilled water or dry. Prolonged storage of the probe in the wet state is not suitable due to the possibility of microbial contamination of the storage water.

Sensor cleaning

Hot water is usually sufficient to remove salts from the probe cell. The probe can be cleaned with a mild detergent solution and / or dilute nitric acid (1% bw). Immerse the probe in the solution and mix for 2 to 3 minutes. Dilute hydrochloric acid (HCl) or sulfuric acid (H₂SO₄) can also be used for cleaning. If more efficient cleaning is required, use a solution containing 50% isopropanol with HCl. Organic compounds can be easily removed with acetone. Do not rub the platinum electrodes (wires, plates) in the cell with abrasive objects, or otherwise clean them mechanically. They may be damaged and the probe may malfunction! Finally, rinse the cell several times with distilled or deionized water and calibrate.

Table: Sensor specifications

Sensor Name	Cell Constant C	Min. Insertion Depth (mm)	Range
CHS CondiGo 30W	1	30	10 µS/cm to 20 mS/cm
CHS CondiLab 22P	1	22	10 µS/cm to 20 mS/cm
CHS CondiLab 50R	10	50	100 µS/cm to 500 mS/cm
CHS CondiLab 40L	0,1	40	0,1 µS/cm to 1 mS/cm

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The warranty period for sensors is 3 months and does not apply to normal wear and life, which depends on application.