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USER MANUAL

LABORATORY SENSORS FOR pH MEASUREMENT

Sensors (electrodes) overview

This instruction is intended for use with following sensors: *XS pH sensors (supplied separately or with XS instruments), XS testers pH1, pH5, PCS, pX4, pH sensors supplied by Chromservis*

The first use (sensor conditioning)

Before the first use, immerse the pH sensor into pH 4 or pH 7 buffer solution for at least 20 minutes. Never leave the pH sensor in distilled or deionized water. After the conditioning, rinse the sensor with distilled or deionized water. Now is the sensor ready for calibration or pH measurement.

Measurement

Calibrate the pH sensor before measurement of the samples (see below). Rinse the sensor between the samples by distilled or deionized water to avoid cross contamination. **DO NOT WIPE the sensor** – generated static electricity could cause incorrect measurements. Remove excess liquid using pure filtration paper. **Electrodes with plastic shaft never do not insert into organic solvents/matrixes.**

Sensor calibration

Sensors have to be periodically calibrated by means of two or three pH buffers. The interval of calibrations is dependent on required precision and accuracy (usually once a day, week or month). For example, if you measure values below than pH 7, calibrate the sensor using pH 4 and pH 7 buffers; for measurements of samples with pH value above 7 calibrate the sensor using pH 7 and pH 10 buffers. After the calibration, make a final check by measuring buffer of pH 7. If the reading is not correct, repeat the calibration. If the calibration buffers were store in the refrigerator, leave them at room temperature to equilibrate (correct values are at the buffer label). The pH meter check should be done before each sample series measurement.

Storage

Sensor must be stored wet, preferably in storage solution supplied with the sensor. This will ensure the longest possible sensor lifetime. Other option is to store the sensor in storage solution (CHS-3KCL-050). Rinse the sensor before the measurement. If there are white crystals, carefully remove them. Remove excess liquid using pure filtration paper. **Never store the pH sensor in distilled and deionized water, protein solutions or solution with fluorides. Do not let it dry. The sensors get older even though they are not used.**

Filling the sensor (not applicable at maintenance free sensors)

Fill the sensor with recommended electrolyte through filling hole. Before you take a measurement, ensure that the sensor includes enough electrolyte and the hole is not covered or capped. Proper electrolyte level must ensure that diaphragm and Ag/AgCl wire is below it. Only use electrolytes recommended by the manufacturer (they may differ).

Sensor cleaning

Grease and fat removal

Rinse the sensor with weak solution of detergent or methanol for 5 minutes. Rinse the sensor with distilled or deionized water before use.

Bacteria and organic deposits removal

Immerse the sensor in diluted solution (1:10) of sodium hypochlorite. Rinse the sensor with distilled or deionized water before use.

Protein deposits removal

Immerse the sensor in Protein Cleaning Solution (CHS-PCS-050) for 5 minutes. Rinse the sensor with distilled or deionized water before use.

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The warranty period for sensors is 3 months and does not apply to normal wear and life, which depends on application.

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