

MSPE drawings




Table No. 1: Outer tube

Dimension described	Dimension [mm]	Note
Tube OD	10.3	
Rim OD	12.5	
Rim height	1.8	
Tube total height	41.2	
Tube height (cylinder only)	37.8	

Table No. 2: Inner tube

Dimension described	Dimension [mm]	Note
Rim OD	12.3	
Under rim cylinder OD	11.5	
Tube OD below rim cylinder	8.9	
Tube OD at the bottom	8.4	Frit level
Tube total height	30.0	
Tube height to the bottom	25.6	Frit level
Tube height to luer tip	27.0	

